

Anti-Interferon alpha 2 [2A7] Monoclonal Antibody

AGMB04684

Description

This Anti-Interferon alpha 2 [2A7] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04684
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	IFNA2, IFNA, IFNα2, Interferon alpha-2, IFN-alpha-2, Interferon alpha-A, LeIF A, IFNA2B, IFN-alphaA
Applications:	ELISA
Research Area:	Immunology
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	2A7
Clonality:	Monoclonal Antibody
SwissProt ID:	P01563
Immunogen:	Recombinant human IFN α 2a
Gene ID:	3440

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: IFNA2

Host Species: Mouse

Isotype: IgG1

Purification: Affinity Purified

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 0.09% sodium azide as a preservative.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	ELISA	1:10000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.