

Anti-NAV1 Rat [7F12] Monoclonal Antibody

AGMB04706

Description

This Anti-NAV1 Rat [7F12] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04706
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	NAV1, Pore membrane and/or filament-interacting-like protein 3, Steerin-1, Unc-53 homolog 1, unc53H1
Applications:	WB ICC/IF
Research Area:	-
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	7F12
Clonality:	Monoclonal Antibody
SwissProt ID:	Q8NEY1
Immunogen:	Recombinant NAV1
Gene ID:	89796

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name:	NAV1
Host Species:	Rat
Isotype:	IgG2a
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	Calculated MW: 202 kDa, Observed MW: 202 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 0.09% sodium azide as a preservative.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	IF	1:50-1:200
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.	