

Anti-FCAMR [10C4] Monoclonal Antibody

AGMB04707

Description

This Anti-FCAMR [10C4] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04707
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	FCA/MR, FKSG87, Fc alpha/mu receptor, CD351
Applications:	IHC-F IHC-P
Research Area:	-
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	10C4
Clonality:	Monoclonal Antibody
SwissProt ID:	Q8WWV6
Immunogen:	Recombinant human Fc γ 2bR
Gene ID:	83953
Gene Name:	FCAMR

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Host Species:	Mouse
Isotype:	IgG2a
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	Calculated MW: 57 kDa, Observed MW: 57 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 0.09% sodium azide as a preservative.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	IHC	1:500

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.