

# Anti-KLH [9G11] Monoclonal Antibody

## AGMB04712

### Description

---

This Anti-KLH [9G11] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

|                       |                                               |
|-----------------------|-----------------------------------------------|
| <b>SKU:</b>           | AGMB04712                                     |
| <b>Contents:</b>      | 50µl, 100µl<br>Bradford Reagent: 1 vial (2ml) |
| <b>Synonyms:</b>      | Keyhole limpet hemocyanin                     |
| <b>Applications:</b>  | <b>ELISA</b>                                  |
| <b>Research Area:</b> | Epigenetics                                   |
| <b>Form:</b>          | Liquid                                        |

### Antibody Data

---

|                      |                           |
|----------------------|---------------------------|
| <b>Reactivity:</b>   | Isotype Control           |
| <b>Clone:</b>        | 9G11                      |
| <b>Clonality:</b>    | Monoclonal Antibody       |
| <b>SwissProt ID:</b> | Q10583                    |
| <b>Immunogen:</b>    | keyhole limpet hemocyanin |
| <b>Gene ID:</b>      | -                         |
| <b>Gene Name:</b>    | -                         |

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Host Species:** Mouse

**Isotype:** IgG2b

**Purification:** Affinity Purified

**Conjugated:** Unconjugated

**Modification:** Unmodified

**Molecular Weight:** -

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in PBS containing 0.09% sodium azide as a preservative.

| Antibody Dilution Ratio: | Application | Antibody Dilution Ratio |
|--------------------------|-------------|-------------------------|
|                          | ELISA       | 1:10000                 |

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.