

Anti-CD74 [R09-8J3] Monoclonal Antibody

AGMB04713

Description

This Anti-CD74 [R09-8J3] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04713
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	CD74, DHLA, HLA class II histocompatibility antigen gamma chain, HLA-DR antigens-associated invariant chain, Ia antigen-associated invariant chain, Ii, p33, CD74
Applications:	WB IHC-P ICC/IF FC
Research Area:	Immunology
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R09-8J3
Clonality:	Monoclonal Antibody
SwissProt ID:	P04233
Immunogen:	A synthesized peptide derived from human CD74
Gene ID:	972

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: CD74

Host Species: Rabbit

Isotype: IgG

Purification: Affinity Chromatography

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW: 34 kDa, Observed MW: 34 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:2000
	IHC	1:100-1:200
	IF	1:50-1:200
	FC	1:50-1:100

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.