

Anti-ERBB2 (D769H Mutant) [7Z5-K7-U6] Monoclonal Antibody

AGMB04764

Description

This Anti-ERBB2 (D769H Mutant) [7Z5-K7-U6] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04764
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	ERBB2,HER2,Neu,c-Neu
Applications:	WB ICC/IF ELISA
Research Area:	Tags & Cell Markers
Form:	Liquid

Antibody Data

Reactivity:	Vertebrates
Clone:	7Z5-K7-U6
Clonality:	Monoclonal Antibody
SwissProt ID:	-
Immunogen:	A synthetic peptide from the internal region of ERBB2 which includes the mutation of D769H, human origin.
Gene ID:	-

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: ERBB2,HER2,Neu,c-Neu

Host Species: Mouse

Isotype: IgG

Purification: Purified from ascites

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: -

Preparation & Storage

Storage: Store at -20°C. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS (without Mg²⁺ and Ca²⁺), pH 7.4, 150 mM NaCl, 50% glycerol

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	ELISA	1:1000-1:2000
	WB	1:500-1:1000
	ICC/IF	1:25-1:50

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.