

# Anti-PIK3CA (Q546K Mutant) [4F4-P3-Y9] Monoclonal Antibody

## AGMB04811

### Description

---

This Anti-PIK3CA (Q546K Mutant) [4F4-P3-Y9] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

|                       |                                               |
|-----------------------|-----------------------------------------------|
| <b>SKU:</b>           | AGMB04811                                     |
| <b>Contents:</b>      | 50µl, 100µl<br>Bradford Reagent: 1 vial (2ml) |
| <b>Synonyms:</b>      | PIK3CA                                        |
| <b>Applications:</b>  | <b>WB</b> <b>IHC-P</b> <b>ELISA</b>           |
| <b>Research Area:</b> | Signal Transduction                           |
| <b>Form:</b>          | Liquid                                        |

### Antibody Data

---

|                      |                                                                                |
|----------------------|--------------------------------------------------------------------------------|
| <b>Reactivity:</b>   | Vertebrates                                                                    |
| <b>Clone:</b>        | 4F4-P3-Y9                                                                      |
| <b>Clonality:</b>    | Monoclonal Antibody                                                            |
| <b>SwissProt ID:</b> | -                                                                              |
| <b>Immunogen:</b>    | A synthetic peptide from the internal region of PIK3CA(Q546K) of human origin. |
| <b>Gene ID:</b>      | -                                                                              |

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Gene Name:** PIK3CA

**Host Species:** Mouse

**Isotype:** -

**Purification:** Purified from ascites

**Conjugated:** Unconjugated

**Modification:** Unmodified

**Molecular Weight:** -

## Preparation & Storage

---

**Storage:** Store at -20°C. Avoid freeze/thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150 mM NaCl, 50% glycerol

| Antibody Dilution Ratio: | Application | Antibody Dilution Ratio |
|--------------------------|-------------|-------------------------|
|                          | ELISA       | 1:1000-1:5000           |
|                          | WB          | 1:100-1:1000            |
|                          | IHC-P       | 1:50-1:100              |

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.