

Anti-PIK3CA (D549N Mutant) [9M2-R3-S8] Monoclonal Antibody

AGMB04816

Description

This Anti-PIK3CA (D549N Mutant) [9M2-R3-S8] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04816
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	PIK3CA, CLOVE, CWS5, MCAP, MCM, MCMTc, PI3K, p110-α
Applications:	WB IHC-P ELISA
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Vertebrates
Clone:	9M2-R3-S8
Clonality:	Monoclonal Antibody
SwissProt ID:	-
Immunogen:	A synthetic peptide from the internal region of PIK3CA which includes the mutation of D549N, human origin.
Gene ID:	-

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: PIK3CA, CLOVE, CWS5, MCAP, MCM, MCMTC, PI3K, p110-alpha

Host Species: Mouse

Isotype: IgG

Purification: Purified from ascites

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: -

Preparation & Storage

Storage: Store at -20°C. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS (without Mg²⁺ and Ca²⁺), pH 7.4, 150 mM NaCl, 50% glycerol

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	IHC-P	1:50-1:100
	ELISA	1:1000-1:5000
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.	