

# Anti-Ras (A146V Mutant) [1X6-S2-T2] Monoclonal Antibody

## AGMB04835

### Description

---

This Anti-Ras (A146V Mutant) [1X6-S2-T2] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

|                       |                                               |
|-----------------------|-----------------------------------------------|
| <b>SKU:</b>           | AGMB04835                                     |
| <b>Contents:</b>      | 50µl, 100µl<br>Bradford Reagent: 1 vial (2ml) |
| <b>Synonyms:</b>      | H-Ras, K-Ras, N-Ras                           |
| <b>Applications:</b>  | <b>WB</b> <b>IHC-P</b> <b>ELISA</b>           |
| <b>Research Area:</b> | Cell Biology                                  |
| <b>Form:</b>          | Liquid                                        |

### Antibody Data

---

|                      |                                                                                                         |
|----------------------|---------------------------------------------------------------------------------------------------------|
| <b>Reactivity:</b>   | Vertebrates                                                                                             |
| <b>Clone:</b>        | 1X6-S2-T2                                                                                               |
| <b>Clonality:</b>    | Monoclonal Antibody                                                                                     |
| <b>SwissProt ID:</b> | -                                                                                                       |
| <b>Immunogen:</b>    | A synthetic peptide from the internal region of Ras which includes the mutation of A146V, human origin. |
| <b>Gene ID:</b>      | -                                                                                                       |

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Gene Name:** H-Ras, K-Ras, N-Ras

**Host Species:** Mouse

**Isotype:** IgG

**Purification:** Purified from ascites

**Conjugated:** Unconjugated

**Modification:** Unmodified

**Molecular Weight:** -

## Preparation & Storage

---

**Storage:** Store at -20°C. Avoid freeze/thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150 mM NaCl, 50% glycerol

| Antibody Dilution Ratio:                  | Application                                                                                                                                                                                                                                             | Antibody Dilution Ratio |
|-------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------|
|                                           | WB                                                                                                                                                                                                                                                      | 1:500-1:1000            |
|                                           | IHC-P                                                                                                                                                                                                                                                   | 1:50-1:100              |
|                                           | ELISA                                                                                                                                                                                                                                                   | 1:1000-1:5000           |
| <b>Protein Quantification (Optional):</b> | To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <a href="https://www.assaygenie.com/bradford-protein-assay-protocol/">https://www.assaygenie.com/bradford-protein-assay-protocol/</a> to view the full protocol. |                         |