

Anti-Tau [3M8-G1-K1] Monoclonal Antibody

AGMB04853

Description

This Anti-Tau [3M8-G1-K1] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04853
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	TAU, MSTD, PPND, DDPAC, MAPTL, MTBT1, MTBT2, tau-40, FTDP-17, PPP1R103, Tau-PHF6, Tau
Applications:	
Research Area:	Neuroscience
Form:	Liquid

Antibody Data

Reactivity:	Mouse, Rat
Clone:	3M8-G1-K1
Clonality:	Monoclonal Antibody
SwissProt ID:	P10636
Immunogen:	Peptide
Gene ID:	4137

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: MAPT
Host Species: Mouse
Isotype: IgG1
Purification: Protein A
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: Calculated MW: 80 kDa, Observed MW:50-65 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS, Glycerol and BSA

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.