

Anti-MCM7/PRL [1H6-B7-M2] Monoclonal Antibody

AGMB04855

Description

This Anti-MCM7/PRL [1H6-B7-M2] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04855
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	MCM7, CDC47, MCM2, DNA replication licensing factor MCM7, CDC47 homolog, P1.1-MCM3, CDABP0042, CDC 47, CDC47,
Applications:	
Research Area:	Epigenetics
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	1H6-B7-M2
Clonality:	Monoclonal Antibody
SwissProt ID:	P33993
Immunogen:	Peptide
Gene ID:	4176

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: MCM7
Host Species: Mouse
Isotype: IgG1
Purification: Protein A
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: Calculated MW: 81 kDa, Observed MW:85 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS, Glycerol and BSA

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.