

# Anti-CD81 [8M2-M8-Z3] Monoclonal Antibody

## AGMB04858

### Description

---

This Anti-CD81 [8M2-M8-Z3] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

<b>SKU:</b>	AGMB04858
<b>Contents:</b>	50µl, 100µl Bradford Reagent: 1 vial (2ml)
<b>Synonyms:</b>	26 kDa cell surface protein TAPA 1, 26 kDa cell surface protein TAPA-1, 26 kDa cell surface protein TAPA1,
<b>Applications:</b>	<b>ICC/IF</b> <b>FC</b>
<b>Research Area:</b>	Immunology
<b>Form:</b>	Liquid

### Antibody Data

---

<b>Reactivity:</b>	Human
<b>Clone:</b>	8M2-M8-Z3
<b>Clonality:</b>	Monoclonal Antibody
<b>SwissProt ID:</b>	P60033
<b>Immunogen:</b>	Peptide
<b>Gene ID:</b>	975

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Gene Name:** CD81  
**Host Species:** Mouse  
**Isotype:** IgG1  
**Purification:** Protein A  
**Conjugated:** Unconjugated  
**Modification:** Unmodified  
**Molecular Weight:** -

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
 Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in PBS, Glycerol and BSA

<b>Antibody Dilution Ratio:</b>	<b>Application</b>	<b>Antibody Dilution Ratio</b>
	ICC/IF	1:50
	FC	1:50-1:100

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.