

Anti-mTOR [9E2-B3-V8] Monoclonal Antibody

AGMB04861

Description

This Anti-mTOR [9E2-B3-V8] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04861
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	MTOR, FRAP, FRAP1, FRAP2, RAFT1, RAPT1, Serine/threonine-protein kinase mTOR, FK506-binding protein 12-rapamycin complex-associated protein 1, FKBP12-rapamycin complex-associated protein, Mammalian target of rapamycin, mTOR, Mechanistic tar
Applications:	WB
Research Area:	Cell Biology
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	9E2-B3-V8
Clonality:	Monoclonal Antibody
SwissProt ID:	P42345
Immunogen:	Peptide

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene ID: 2475

Gene Name: MTOR

Host Species: Mouse

Isotype: IgG1

Purification: Protein A

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW: 289 kDa, Observed MW:289 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS, Glycerol and BSA

	Application	Antibody Dilution Ratio
Antibody Dilution Ratio:	WB	1:500-1:1000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.