

Anti-Caveolin-1 [2F4-U4-G7] Monoclonal Antibody

AGMB04863

Description

This Anti-Caveolin-1 [2F4-U4-G7] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04863
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	CAV1, CAV, Caveolin-1, BSCL3, CAV, CAV1,
Applications:	WB
Research Area:	Cardiovascular
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse
Clone:	2F4-U4-G7
Clonality:	Monoclonal Antibody
SwissProt ID:	Q03135
Immunogen:	Peptide
Gene ID:	857
Gene Name:	CAV1

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Host Species: Mouse

Isotype: IgG1

Purification: Protein G

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW: 20 kDa, Observed MW:20 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS, Glycerol and BSA

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.