

Anti-HIP1 [1L4-D1-V4] Monoclonal Antibody

AGMB04872

Description

This Anti-HIP1 [1L4-D1-V4] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04872
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	2610109B09Rik, A930014B11Rik, E130315I21Rik,
Applications:	WB ICC/IF
Research Area:	Neuroscience
Form:	Liquid

Antibody Data

Reactivity:	Mouse
Clone:	1L4-D1-V4
Clonality:	Monoclonal Antibody
SwissProt ID:	O00291
Immunogen:	Peptide
Gene ID:	3092
Gene Name:	HIP1

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Host Species: Mouse

Isotype: IgG1

Purification: Protein A

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW: 116 kDa, Observed MW:116 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS, Glycerol and BSA

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	ICC/IF	1:50

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.