

Anti-Pleckstrin [8E8-C3-P3] Monoclonal Antibody

AGMB04886

Description

This Anti-Pleckstrin [8E8-C3-P3] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04886
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	PLEK, P47, Platelet 47 kDa protein
Applications:	WB IHC-P ICC/IF
Research Area:	Cardiovascular
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	8E8-C3-P3
Clonality:	Monoclonal Antibody
SwissProt ID:	P08567
Immunogen:	Peptide
Gene ID:	5341
Gene Name:	PLEK

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Host Species: Mouse

Isotype: IgG1

Purification: Protein G

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW: 40 kDa, Observed MW:47 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS, Glycerol and BSA

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	IHC-P	1:50-1:100
	ICC/IF	1:50

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.