

Anti-GRP94 [6F6-W3-U2] Monoclonal Antibody

AGMB04890

Description

This Anti-GRP94 [6F6-W3-U2] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04890
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	ECGP, GP96, TRA1, GRP94, HEL35, HEL-S-125m, Grp94
Applications:	WB
Research Area:	Tags & Cell Markers
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	6F6-W3-U2
Clonality:	Monoclonal Antibody
SwissProt ID:	P14625
Immunogen:	Peptide
Gene ID:	7184
Gene Name:	HSP90B1

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Host Species: Mouse

Isotype: IgG1

Purification: Protein G

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW: 95 kDa, Observed MW:100 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS, Glycerol and BSA

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.