

# Anti-Parvalbumin [1D9-X6-C1] Monoclonal Antibody

AGMB04892

## Description

---

This Anti-Parvalbumin [1D9-X6-C1] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

## Product Information

---

<b>SKU:</b>	AGMB04892
<b>Contents:</b>	50µl, 100µl Bradford Reagent: 1 vial (2ml)
<b>Synonyms:</b>	Parvalbumin alpha, PVALB, D22S749, MGC116759, Parvalbumin alpha,
<b>Applications:</b>	<b>WB</b> <b>IHC-P</b>
<b>Research Area:</b>	Neuroscience
<b>Form:</b>	Liquid

## Antibody Data

---

<b>Reactivity:</b>	Mouse, Rat
<b>Clone:</b>	1D9-X6-C1
<b>Clonality:</b>	Monoclonal Antibody
<b>SwissProt ID:</b>	P20472
<b>Immunogen:</b>	Peptide
<b>Gene ID:</b>	5816
<b>Gene Name:</b>	PVALB

### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Host Species:** Mouse

**Isotype:** IgG1

**Purification:** Protein G

**Conjugated:** Unconjugated

**Modification:** Unmodified

**Molecular Weight:** Calculated MW: 14 kDa, Observed MW:14 kDa

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in PBS, Glycerol and BSA

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	IHC-P	1:50-1:100

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.