

Anti-CAPG/MCP [3N9-A7-T5] Monoclonal Antibody

AGMB04896

Description

This Anti-CAPG/MCP [3N9-A7-T5] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04896
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	AFCP, Capg, gCap39, HEL S 66, mbh1, MCP, Actin capping protein GCAP39, Actin regulatory protein CAP G, Actin regulatory protein CAP-G,
Applications:	WB ICC/IF
Research Area:	Immunology
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	3N9-A7-T5
Clonality:	Monoclonal Antibody
SwissProt ID:	P40121
Immunogen:	Peptide
Gene ID:	822

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: CAPG
Host Species: Mouse
Isotype: IgG1
Purification: Protein G
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: Calculated MW: 40 kDa, Observed MW:40 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS, Glycerol and BSA

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	ICC/IF	1:50

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.