

Anti-Mucin 5AC [4X3-S2-K6] Monoclonal Antibody

AGMB04909

Description

This Anti-Mucin 5AC [4X3-S2-K6] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04909
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	MUC5AC, Apomucin, Gastric mucin, LeB, Mucin 5AC oligomeric mucus/gel forming, TBM
Applications:	IHC-P
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	4X3-S2-K6
Clonality:	Monoclonal Antibody
SwissProt ID:	P98088
Immunogen:	Peptide
Gene ID:	4586

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: MUC5AC
Host Species: Mouse
Isotype: IgG1
Purification: Protein A
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS, Glycerol and BSA

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	IHC-P	1:50-1:100

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.