

Anti-RNA Helicase A [9I6-W3-V4] Monoclonal Antibody

AGMB04911

Description

This Anti-RNA Helicase A [9I6-W3-V4] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04911
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	DHX9/RNA Helicase A, LKP, RHA, DDX9, NDH2, NDHII, ATP dependent RNA helicase A, ATP-dependent RNA helicase A, DDX 9,
Applications:	WB IHC-P ICC/IF
Research Area:	Epigenetics
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	9I6-W3-V4
Clonality:	Monoclonal Antibody
SwissProt ID:	Q08211
Immunogen:	Peptide
Gene ID:	1660

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: DHX9
Host Species: Mouse
Isotype: IgG1
Purification: Protein A
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: Calculated MW: 141 kDa, Observed MW:141kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS, Glycerol and BSA

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	IHC-P	1:50-1:100
	ICC/IF	1:50
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.	