

Anti-PDCD4 [9T7-M7-I9] Monoclonal Antibody

AGMB04919

Description

This Anti-PDCD4 [9T7-M7-I9] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04919
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	Death up-regulated gene protein, Dug, H731, Ma3, Neoplastic transformation inhibitor, Nuclear antigen H731, Pdcd4, Programmed cell death protein 4, Tis
Applications:	WB IHC-P ICC/IF
Research Area:	Cell Biology
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	9T7-M7-I9
Clonality:	Monoclonal Antibody
SwissProt ID:	Q53EL6
Immunogen:	Peptide
Gene ID:	27250

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: PDCD4

Host Species: Mouse

Isotype: IgG1

Purification: Protein A

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW: 52 kDa, Observed MW:56 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS, Glycerol and BSA

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	IHC-P	1:50-1:100
	ICC/IF	1:50
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.	