

Anti-SLAMF7/CS1 [8H9-Z3-H8] Monoclonal Antibody

AGMB04933

Description

This Anti-SLAMF7/CS1 [8H9-Z3-H8] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04933
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	CD319, 19A, CS1, SLAMF7, SLAM family member 7
Applications:	ICC/IF FC
Research Area:	Immunology
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	8H9-Z3-H8
Clonality:	Monoclonal Antibody
SwissProt ID:	Q9NQ25
Immunogen:	Peptide
Gene ID:	57823
Gene Name:	SLAMF7

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Host Species: Mouse
Isotype: IgG2b
Purification: Protein G
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS, Glycerol and BSA

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	ICC/IF	1:50
	FC	1:50-1:100

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.