

Anti-Pepsinogen II/PGC [1M9-T5-G3] Monoclonal Antibody

AGMB04958

Description

This Anti-Pepsinogen II/PGC [1M9-T5-G3] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04958
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	PGC, Gastricsin, Pepsinogen C
Applications:	WB IHC-F IHC-P ICC/IF
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	1M9-T5-G3
Clonality:	Monoclonal Antibody
SwissProt ID:	P20142
Immunogen:	Recombinant human PGC protein
Gene ID:	5225
Gene Name:	PGC

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Host Species: Mouse

Isotype: IgG2b

Purification: Affinity Purified by Protein A

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW: 42kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 0.01M PBS ,PH7.4, with 0.02% Proclin300.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:2000
	IHC	1:100-1:500
	IF	1:100-1:500
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.	