

HRP (Mouse)& AP (Rabbit) IHC Detection Kit (Dual Staining)

AGMB04967

Description

This HRP (Mouse)& AP (Rabbit) IHC Detection Kit (Dual Staining) is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04967
Contents:	100-200T, 300-600T Bradford Reagent: 1 vial (2ml)
Synonyms:	-
Applications:	IHC-P
Research Area:	Immunology
Form:	Liquid

Antibody Data

Reactivity:	Mouse, Rabbit
Clone:	-
Clonality:	Polyclonal Antibody
SwissProt ID:	-
Immunogen:	-
Gene ID:	-
Gene Name:	-

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Host Species: -
Isotype: -
Purification: -
Conjugated: HRP/AP
Modification: -
Molecular Weight: -

Preparation & Storage

Storage: Store at 2~8°C. Avoid freezing. Valid for 12 months. The reagent is ready-to-use. After each use, it should be put back into the refrigerator at 2~8°C
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer:

Antibody Dilution Ratio: -

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.