

GenieSignal Multiplex IHC Detection Kit (Sextuple)

AGMB04973

Description

This GenieSignal Multiplex IHC Detection Kit (Sextuple) is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04973
Contents:	20T, 50T, 100T Bradford Reagent: 1 vial (2ml)
Synonyms:	-
Applications:	IHC-P mIHC
Research Area:	Immunology
Form:	Liquid

Antibody Data

Reactivity:	Mouse, Rabbit
Clone:	-
Clonality:	Polyclonal Antibody
SwissProt ID:	-
Immunogen:	-
Gene ID:	-
Gene Name:	-

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Host Species: Goat
Isotype: -
Purification: -
Conjugated: HRP
Modification: -
Molecular Weight: -

Preparation & Storage

Storage: Please refer to the Instructions for details
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Reagent 1: Goat anti-Mouse/Rabbit HRP polymer Reagent 2: TSA Enhancer
 Reagent 3: TSA-480 Dye (Ex: 450nm / Em: 480nm) Reagent 4: TSA-520 Dye (Ex: 488nm / Em: 520nm)
 Reagent 5: TSA-570 Dye (Ex: 555nm / Em: 570nm) Reagent 6: TSA-620 Dye (Ex: 590nm / Em: 620nm)
 Reagent 7: TSA-690 Dye (Ex: 630nm / Em: 690nm) Reagent 8: TSA-780 Dye (Ex: 750nm / Em: 780nm)

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	TSA Enhancer : TSA Dye	1 : 500

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.