

GenieSignal Multiplex IHC Detection Kit (For Frozen section /Cell slides , Triple)

AGMB04975

Description

This GenieSignal Multiplex IHC Detection Kit (For Frozen section /Cell slides, Triple) is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04975
Contents:	20T, 50T, 100T Bradford Reagent: 1 vial (2ml)
Synonyms:	-
Applications:	IHC-P IHC-F ICC/IF mIHC
Research Area:	Immunology
Form:	Liquid

Antibody Data

Reactivity:	Mouse, Rabbit
Clone:	-
Clonality:	Polyclonal Antibody
SwissProt ID:	-
Immunogen:	-
Gene ID:	-

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: -
Host Species: Goat
Isotype: -
Purification: -
Conjugated: HRP
Modification: -
Molecular Weight: -

Preparation & Storage

Storage: Store at -20°C .Avoid freeze/thaw cycles.Avoid light.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Reagent1:TSA-520 Dye(Ex: 488nm / Em: 520nm) Reagent2:TSA-570 Dye(Ex: 555nm / Em: 570nm) Reagent3:TSA-690 Dye(Ex: 630nm / Em: 690nm)
 Reagent4:TSA Enhancer Reagent5:Antibody Eluent Reagent6:Goat anti-Mouse/Rabbit HRP polymer

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	TSA Enhancer : TSA Dye	1 : 500

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.