

# Sensitive ECL Western Blotting Substrate [Medium-pg Level]

AGMB04976

## Description

---

This Sensitive ECL Western Blotting Substrate [Medium-pg Level] is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

## Product Information

---

<b>SKU:</b>	AGMB04976
<b>Contents:</b>	50mL, 100mL, 500mL Bradford Reagent: 1 vial (2ml)
<b>Synonyms:</b>	-
<b>Applications:</b>	<b>WB</b>
<b>Research Area:</b>	Immunology
<b>Form:</b>	Liquid

## Antibody Data

---

<b>Reactivity:</b>	Species-Independent
<b>Clone:</b>	-
<b>Clonality:</b>	-
<b>SwissProt ID:</b>	-
<b>Immunogen:</b>	-
<b>Gene ID:</b>	-
<b>Gene Name:</b>	-

### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Host Species:** -

**Isotype:** -

**Purification:** -

**Conjugated:** -

**Modification:** -

**Molecular Weight:** -

## Preparation & Storage

---

**Storage:** Store at 2-8°C..Avoid light.Avoid freezing.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** -

**Antibody Dilution Ratio:** -

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.