

Hypersensitive ECL Western Blotting Substrate [Low-pg Level]

AGMB04977

Description

This Hypersensitive ECL Western Blotting Substrate [Low-pg Level] is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04977
Contents:	50mL, 100mL, 500mL Bradford Reagent: 1 vial (2ml)
Synonyms:	-
Applications:	WB
Research Area:	Immunology
Form:	Liquid

Antibody Data

Reactivity:	Species-Independent
Clone:	-
Clonality:	-
SwissProt ID:	-
Immunogen:	-
Gene ID:	-
Gene Name:	-

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Host Species: -

Isotype: -

Purification: -

Conjugated: -

Modification: -

Molecular Weight: -

Preparation & Storage

Storage: Store at 2-8°C..Avoid light.Avoid freezing.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: -

Antibody Dilution Ratio: -

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.