

Anti-Acetyl-Alpha Tubulin (Lys40) (6B5) [6B5-1G7-3F9] Monoclonal Antibody

AGMB04986

Description

This Anti-Acetyl-Alpha Tubulin (Lys40) (6B5) [6B5-1G7-3F9] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04986
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	Alpha Tubulin 2, B Alpha 1, H2 Alpha, Hum a tub 1, TUBA1, TUBA1A, Tubulin Alpha 1Chain
Applications:	WB IHC-P
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	6B5-1G7-3F9
Clonality:	Monoclonal Antibody
SwissProt ID:	P68363
Immunogen:	Synthetic peptide conjugated to KLH.
Gene ID:	10376

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: TUBA1B
Host Species: Mouse
Isotype: IgG1
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Acetylated
Molecular Weight: Calculated MW: 50 kDa, Observed MW: 50 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	IHC	1:50-1:100

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.