

Anti-Acetyl-Histone H3 (Lys4) [9S1-O2-A4] Monoclonal Antibody

AGMB05003

Description

This Anti-Acetyl-Histone H3 (Lys4) [9S1-O2-A4] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05003
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	H3K4ac, H3/j, H3C1, H3C2, H3C3, H3C4, H3C6, H3C7, H3C8, H3FJ, H3C10, H3C11, HIST1H3J
Applications:	
Research Area:	Epigenetics
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	9S1-O2-A4
Clonality:	Monoclonal Antibody
SwissProt ID:	P68431
Immunogen:	Peptide
Gene ID:	8350

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: H3C1

Host Species: Mouse

Isotype: IgG

Purification: Protein G

Conjugated: Unconjugated

Modification: Acetylated

Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS, Glycerol and BSA

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.