

# Anti-Acetyl-Histone H2A (Lys15) [6V5-Z6-T8] Monoclonal Antibody

## AGMB05008

### Description

---

This Anti-Acetyl-Histone H2A (Lys15) [6V5-Z6-T8] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

<b>SKU:</b>	AGMB05008
<b>Contents:</b>	50µl, 100µl Bradford Reagent: 1 vial (2ml)
<b>Synonyms:</b>	H2AK15ac, Histone H2A.2 ; Histone H2A/a ; Histone H2A/m
<b>Applications:</b>	<b>WB</b>
<b>Research Area:</b>	Epigenetics
<b>Form:</b>	Liquid

### Antibody Data

---

<b>Reactivity:</b>	Mouse, Rat
<b>Clone:</b>	6V5-Z6-T8
<b>Clonality:</b>	Monoclonal Antibody
<b>SwissProt ID:</b>	P04908
<b>Immunogen:</b>	Peptide
<b>Gene ID:</b>	3012
<b>Gene Name:</b>	H2AC4

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Host Species:** Mouse  
**Isotype:** IgG  
**Purification:** Protein G  
**Conjugated:** Unconjugated  
**Modification:** Acetylated  
**Molecular Weight:** -

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
 Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in PBS, Glycerol and BSA

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.