

Anti-Butyryl-Histone H2B (Lys20) [6Q8-G5-B2] Monoclonal Antibody

AGMB05029

Description

This Anti-Butyryl-Histone H2B (Lys20) [6Q8-G5-B2] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05029
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	H2BK11ac, Histone H2B.1 A, Histone H2B.a (H2B/a), Histone H2B.g (H2B/g), Histone H2B.h (H2B/h), Histone H2B.k (H2B/k), Histone H2B.l (H2B/l), GL105, H2B GL105, H2B histone family member O,
Applications:	WB
Research Area:	Epigenetics
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	6Q8-G5-B2
Clonality:	Monoclonal Antibody
SwissProt ID:	P62807
Immunogen:	Peptide
Gene ID:	-

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: H2BC4

Host Species: Mouse

Isotype: IgG

Purification: Protein G

Conjugated: Unconjugated

Modification: Butyrylated

Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS, Glycerol and BSA

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.