

# Anti-Phospho-RIP3 (Ser232) [R66-4K-6] Monoclonal Antibody

## AGMB05037

### Description

---

This Anti-Phospho-RIP3 (Ser232) [R66-4K-6] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

<b>SKU:</b>	AGMB05037
<b>Contents:</b>	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
<b>Synonyms:</b>	RIP3, RIP like protein kinase 3, RIP 3, Receptor-interacting serine/threonine-protein kinase 3, Receptor-interacting protein 3
<b>Applications:</b>	<b>WB</b>
<b>Research Area:</b>	Cell Biology
<b>Form:</b>	Liquid

### Antibody Data

---

<b>Reactivity:</b>	Mouse
<b>Clone:</b>	R66-4K-6
<b>Clonality:</b>	Monoclonal Antibody
<b>SwissProt ID:</b>	Q9Y572/Q9QZL0
<b>Immunogen:</b>	A synthesized peptide derived from human Phospho-RIP3 (S232)
<b>Gene ID:</b>	11035/56532

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Gene Name:** RIPK3/Ripk3  
**Host Species:** Rabbit  
**Isotype:** IgG  
**Purification:** Affinity Chromatography  
**Conjugated:** Unconjugated  
**Modification:** Phosphorylated  
**Molecular Weight:** Calculated MW:53/57kDa, Observed MW:50kDa

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
 Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

<b>Antibody Dilution Ratio:</b>	<b>Application</b>	<b>Antibody Dilution Ratio</b>
	WB	1:500-1:1000

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.