

# Anti-Phospho-CDC37 (Ser13) [R04-4I8] Monoclonal Antibody

## AGMB05054

### Description

---

This Anti-Phospho-CDC37 (Ser13) [R04-4I8] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

|                       |                                                                                                     |
|-----------------------|-----------------------------------------------------------------------------------------------------|
| <b>SKU:</b>           | AGMB05054                                                                                           |
| <b>Contents:</b>      | 20µl, 50µl, 100µl<br>Bradford Reagent: 1 vial (2ml)                                                 |
| <b>Synonyms:</b>      | CDC37, CDC37A, Hsp90 co-chaperone Cdc37, Hsp90 chaperone protein kinase-targeting subunit, p50Cdc37 |
| <b>Applications:</b>  |                  |
| <b>Research Area:</b> | Cell Biology                                                                                        |
| <b>Form:</b>          | Liquid                                                                                              |

### Antibody Data

---

|                      |                                                                                       |
|----------------------|---------------------------------------------------------------------------------------|
| <b>Reactivity:</b>   | Human, Mouse, Rat                                                                     |
| <b>Clone:</b>        | R04-4I8                                                                               |
| <b>Clonality:</b>    | Monoclonal Antibody                                                                   |
| <b>SwissProt ID:</b> | Q16543                                                                                |
| <b>Immunogen:</b>    | A synthetic phosphopeptide corresponding to residues surrounding Ser13 of human Cdc37 |
| <b>Gene ID:</b>      | 11140                                                                                 |

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Gene Name:** CDC37

**Host Species:** Rabbit

**Isotype:** IgG

**Purification:** Affinity Purified

**Conjugated:** Unconjugated

**Modification:** Phosphorylated

**Molecular Weight:** Calculated MW: 44 kDa, Observed MW: 44 kDa

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

| Antibody Dilution Ratio: | Application | Antibody Dilution Ratio |
|--------------------------|-------------|-------------------------|
|                          | WB          | 1:1000-1:5000           |

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.