

Anti-Phospho-NRF2 (Ser40) [R03-4A1] Monoclonal Antibody

AGMB05057

Description

This Anti-Phospho-NRF2 (Ser40) [R03-4A1] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05057
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	NFE2L2, NRF2, Nuclear factor erythroid 2-related factor 2, NF-E2-related factor 2, NFE2-related factor 2, HEBP1, Nuclear factor, erythroid derived 2, like 2
Applications:	WB IHC-P ICC/IF FC
Research Area:	Epigenetics
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R03-4A1
Clonality:	Monoclonal Antibody
SwissProt ID:	Q16236
Immunogen:	A synthesized peptide derived from human Phospho-Nrf2 (S40)
Gene ID:	4780

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: NFE2L2
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Chromatography
Conjugated: Unconjugated
Modification: Phosphorylated
Molecular Weight: Calculated MW: 68 kDa, Observed MW: 90 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:5000
	IHC-P	1:100-1:500
	ICC/IF	1:100-1:200
	FC	1:50-1:100

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.