

Anti-Phospho-KIF22 (Ser427) [R06-3E7] Monoclonal Antibody

AGMB05064

Description

This Anti-Phospho-KIF22 (Ser427) [R06-3E7] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05064
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	KIF22, KID, KNSL4, Kinesin-like protein KIF22, Kinesin-like DNA-binding protein, Kinesin-like protein 4
Applications:	
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R06-3E7
Clonality:	Monoclonal Antibody
SwissProt ID:	Q14807
Immunogen:	A synthetic phosphopeptide corresponding to residues surrounding Ser427 of human KIF22
Gene ID:	3835

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: KIF22

Host Species: Rabbit

Isotype: IgG

Purification: Affinity Purified

Conjugated: Unconjugated

Modification: Phosphorylated

Molecular Weight: Calculated MW: 73 kDa, Observed MW: 73 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.