

# Anti-Phospho-GSK3 (alpha + beta) (Tyr279/Tyr216) [R01-5G6] Monoclonal Antibody AGMB05096

## Description

---

This Anti-Phospho-GSK3 (alpha + beta) (Tyr279/Tyr216) [R01-5G6] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

## Product Information

---

<b>SKU:</b>	AGMB05096
<b>Contents:</b>	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
<b>Synonyms:</b>	GSK3α/β, Serine/threonine-protein kinase GSK3A, Serine/threonine-protein kinase GSK3B
<b>Applications:</b>	<b>WB</b> <b>IHC-P</b> <b>ICC/IF</b> <b>IP</b>
<b>Research Area:</b>	Cell Biology
<b>Form:</b>	Liquid

## Antibody Data

---

<b>Reactivity:</b>	Human, Mouse, Rat
<b>Clone:</b>	R01-5G6
<b>Clonality:</b>	Monoclonal Antibody
<b>SwissProt ID:</b>	P49840/P49841
<b>Immunogen:</b>	A synthetic phosphopeptide corresponding to residues surrounding Tyr216 of human GSK3 alpha

### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Gene ID:** 2931/2932

**Gene Name:** GSK3A/GSK3B

**Host Species:** Rabbit

**Isotype:** IgG

**Purification:** Affinity Purified

**Conjugated:** Unconjugated

**Modification:** Phosphorylated

**Molecular Weight:** Calculated MW: 51 kDa, Observed MW: 47-51 kDa

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

**Antibody Dilution Ratio:**

Application	Antibody Dilution Ratio
WB	1:500-1:1000
IHC-P	1:20-1:50
IP	1:20-1:50

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.