

Anti-Phospho-ERK1/2 (Thr202/Tyr204)/ (Thr185/Tyr187) [R04-9E9] Monoclonal Antibody AGMB05130

Description

This Anti-Phospho-ERK1/2 (Thr202/Tyr204)/ (Thr185/Tyr187) [R04-9E9] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05130
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	MAPK3, ERK1, ERT2, ERK-1, PRKM3, P44ERK1, P44MAPK, HS44KDAP, HUMKER1A, p44-ERK1, p44-MAPK, MAPK1, ERK, p38, p40, p41, ERK2, ERT1, ERK-2, MAPK2, PRKM1, PRKM2, P42MAPK, p41mapk, p42-MAPK.
Applications:	WB IP
Research Area:	Neuroscience
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R04-9E9
Clonality:	Monoclonal Antibody
SwissProt ID:	P27361/P28482
Immunogen:	A synthetic phosphopeptide corresponding to residues surrounding Thr202/Tyr204 of human Erk1

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene ID: 5595/5594

Gene Name: MAPK1/MAPK3

Host Species: Rabbit

Isotype: IgG

Purification: Affinity Purified

Conjugated: Unconjugated

Modification: Phosphorylated

Molecular Weight: Calculated MW: 44,42 kDa, Observed MW: 44,42 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:5000
	IP	1:50-1:100

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.