

Anti-Phospho-c-Myc (Thr58) [R09-1M8] Monoclonal Antibody

AGMB05149

Description

This Anti-Phospho-c-Myc (Thr58) [R09-1M8] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05149
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	MYC, BHLHE39, Myc proto-oncogene protein, Class E basic helix-loop-helix protein 39, bHLHe39, Proto-oncogene c-Myc, Transcription factor p64
Applications:	WB ICC/IF FC
Research Area:	Epigenetics
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R09-1M8
Clonality:	Monoclonal Antibody
SwissProt ID:	P01106
Immunogen:	A synthesized peptide derived from human Phospho-c-Myc (T58)
Gene ID:	4609

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: MYC

Host Species: Rabbit

Isotype: IgG

Purification: Affinity Chromatography

Conjugated: Unconjugated

Modification: Phosphorylated

Molecular Weight: Calculated MW: 49 kDa, Observed MW: 57 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:2000
	IF	1:50-1:200
	FC	1:50-1:100

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.