

Anti-Phospho-Presenilin 2/AD5 (Ser330) [R01-5C9] Monoclonal Antibody

AGMB05177

Description

This Anti-Phospho-Presenilin 2/AD5 (Ser330) [R01-5C9] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05177
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	AD3LP, AD5, E5-1, STM-2, AD3L, AD4,
Applications:	WB
Research Area:	Neuroscience
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R01-5C9
Clonality:	Monoclonal Antibody
SwissProt ID:	P49810
Immunogen:	A synthetic phosphopeptide corresponding to residues surrounding Ser330 of human Presenilin 2
Gene ID:	5664

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: PSEN2

Host Species: Rabbit

Isotype: IgG

Purification: Affinity Purified

Conjugated: Unconjugated

Modification: Phosphorylated

Molecular Weight: Calculated MW: 50 kDa, Observed MW: 23 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:5000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.