

Anti-Phospho-beta Arrestin 1 (Ser412) [R05-6F7] Monoclonal Antibody

AGMB05217

Description

This Anti-Phospho-beta Arrestin 1 (Ser412) [R05-6F7] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05217
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	ARB1, ARR1, ARRB1, ARRB1_HUMAN, Arrestin 2, Arrestin beta 1, Arrestin beta-1, Beta-arrestin-1.
Applications:	WB IHC-P
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R05-6F7
Clonality:	Monoclonal Antibody
SwissProt ID:	P49407
Immunogen:	A synthetic phosphopeptide corresponding to residues surrounding Ser412 of human beta Arrestin 1
Gene ID:	408

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: ARRB1

Host Species: Rabbit

Isotype: IgG

Purification: Affinity Purified

Conjugated: Unconjugated

Modification: Phosphorylated

Molecular Weight: Calculated MW: 47 kDa, Observed MW: 50 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:5000
	IHC	1:50-1:100

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.