

# Anti-Phospho-delta 1 Catenin (Thr916) [R07-5A6] Monoclonal Antibody

## AGMB05221

### Description

---

This Anti-Phospho-delta 1 Catenin (Thr916) [R07-5A6] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

|                       |                                                                 |
|-----------------------|-----------------------------------------------------------------|
| <b>SKU:</b>           | AGMB05221                                                       |
| <b>Contents:</b>      | 20µl, 50µl, 100µl<br>Bradford Reagent: 1 vial (2ml)             |
| <b>Synonyms:</b>      | CAS, p120, BCDS2, CTNND, P120CAS, P120CTN, p120(CAS), p120(CTN) |
| <b>Applications:</b>  | <b>WB</b> <b>IHC-P</b> <b>ICC/IF</b>                            |
| <b>Research Area:</b> | Signal Transduction                                             |
| <b>Form:</b>          | Liquid                                                          |

### Antibody Data

---

|                      |                                                                                                      |
|----------------------|------------------------------------------------------------------------------------------------------|
| <b>Reactivity:</b>   | Human                                                                                                |
| <b>Clone:</b>        | R07-5A6                                                                                              |
| <b>Clonality:</b>    | Monoclonal Antibody                                                                                  |
| <b>SwissProt ID:</b> | O60716                                                                                               |
| <b>Immunogen:</b>    | A synthetic phosphopeptide corresponding to residues surrounding Thr916 of human delta 1 Catenin/CAS |
| <b>Gene ID:</b>      | 1500                                                                                                 |

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

|                          |                                              |
|--------------------------|----------------------------------------------|
| <b>Gene Name:</b>        | CTNND1                                       |
| <b>Host Species:</b>     | Rabbit                                       |
| <b>Isotype:</b>          | IgG                                          |
| <b>Purification:</b>     | Affinity Purified                            |
| <b>Conjugated:</b>       | Unconjugated                                 |
| <b>Modification:</b>     | Phosphorylated                               |
| <b>Molecular Weight:</b> | Calculated MW: 108 kDa, Observed MW: 100 kDa |

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

| <b>Antibody Dilution Ratio:</b>           | Application                                                                                                                                                                                                                                             | Antibody Dilution Ratio |
|-------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------|
|                                           | WB                                                                                                                                                                                                                                                      | 1:2000-1:5000           |
|                                           | IHC-P                                                                                                                                                                                                                                                   | 1:100-1:200             |
|                                           | ICC/IF                                                                                                                                                                                                                                                  | 1:50-1:100              |
| <b>Protein Quantification (Optional):</b> | To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <a href="https://www.assaygenie.com/bradford-protein-assay-protocol/">https://www.assaygenie.com/bradford-protein-assay-protocol/</a> to view the full protocol. |                         |