

Anti-Phospho-Crk p38 (Tyr221) [R04-6L8] Monoclonal Antibody

AGMB05261

Description

This Anti-Phospho-Crk p38 (Tyr221) [R04-6L8] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05261
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	CRK, CRKII, p38, Adapter molecule crk, avian sarcoma virus CT10 (v-crk) oncogene homolog, CRK, CRKII, Proto-oncogene c-Crk
Applications:	WB ICC/IF
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R04-6L8
Clonality:	Monoclonal Antibody
SwissProt ID:	P46108
Immunogen:	A synthesized peptide derived from human Phospho-CrkII (Tyr221)
Gene ID:	1398

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: CRK
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Chromatography
Conjugated: Unconjugated
Modification: Phosphorylated
Molecular Weight: Calculated MW: 34 kDa, Observed MW: 42 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

	Application	Antibody Dilution Ratio
Antibody Dilution Ratio:	WB	1:1000-1:2000
	IF	1:50-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.