

Anti-Phospho-LRRK2 (Ser935) [R94-3M-6] Monoclonal Antibody

AGMB05280

Description

This Anti-Phospho-LRRK2 (Ser935) [R94-3M-6] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05280
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	AURA17, Dardarin antibody, , Leucine rich repeat kinase 2, LRRK 2 antibody, LRRK2, LRRK2_HUMAN, PARK 8, PARK8, RIPK7, ROCO 2, ROCO2
Applications:	WB ICC/IF
Research Area:	Neuroscience
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse
Clone:	R94-3M-6
Clonality:	Monoclonal Antibody
SwissProt ID:	Q5S007
Immunogen:	A synthesized peptide derived from human Phospho-LRRK2 (S935)
Gene ID:	120892

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: LRRK2

Host Species: Rabbit

Isotype: IgG

Purification: Affinity Chromatography

Conjugated: Unconjugated

Modification: Phosphorylated

Molecular Weight: Calculated MW:286kDa, Observed MW:286kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:2000
	IF	1:50-1:200
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.	