

Anti-Phospho-RSK1 p90 (Thr359/Ser363) [R33-5H-9] Monoclonal Antibody

AGMB05282

Description

This Anti-Phospho-RSK1 p90 (Thr359/Ser363) [R33-5H-9] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05282
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	RSK, HU-1, MAPKAPK1A, RPS6KA1, p90 RSK1, RPS6K1, RSK 1 p90, S6K alpha 1, 90 kDa ribosomal protein S6 kinase 1 Ribosomal S6 kinase 1,
Applications:	WB IHC-P ICC/IF IP
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R33-5H-9
Clonality:	Monoclonal Antibody
SwissProt ID:	Q15418
Immunogen:	A synthesized peptide derived from human Phospho-RSK1 p90 (T359+S363)

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene ID: 6195

Gene Name: RPS6KA1

Host Species: Rabbit

Isotype: IgG

Purification: Affinity Chromatography

Conjugated: Unconjugated

Modification: Phosphorylated

Molecular Weight: Calculated MW: 83 kDa, Observed MW: 90kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:

Application	Antibody Dilution Ratio
WB	1:1000-1:2000
IHC	1:100-1:200
IF	1:50-1:200
IP	1:50

Protein Quantification (Optional):

To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.