

Anti-Phospho-Tau (Ser214) [R15-2M-5] Monoclonal Antibody

AGMB05306

Description

This Anti-Phospho-Tau (Ser214) [R15-2M-5] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05306
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	TAU, MSTD, PPND, DDPAC, MAPTL, MTBT1, MTBT2, tau-40, FTDP-17, PPP1R103, Tau-PHF6, Tau
Applications:	WB IHC-P
Research Area:	Neuroscience
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse
Clone:	R15-2M-5
Clonality:	Monoclonal Antibody
SwissProt ID:	P10636
Immunogen:	A synthesized peptide derived from human Phospho-Tau (S214)
Gene ID:	4137

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: MAPT
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Chromatography
Conjugated: Unconjugated
Modification: Phosphorylated
Molecular Weight: Calculated MW:79kDa, Observed MW:50-70kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:2000
	IHC	1:100-1:200
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.	