

# Anti-Phospho-HER2 (Tyr1221 /Tyr1222) [R87-5Y-7] Monoclonal Antibody

## AGMB05309

### Description

---

This Anti-Phospho-HER2 (Tyr1221 /Tyr1222) [R87-5Y-7] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

<b>SKU:</b>	AGMB05309
<b>Contents:</b>	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
<b>Synonyms:</b>	NEU, NGL, ERBB2, TKR1, CD340, HER-2, VSCN2, MLN 19, MLN-19, c-ERB2, c-ERB-2, HER-2/neu, p185(erbB2)
<b>Applications:</b>	<b>WB</b>
<b>Research Area:</b>	Tags & Cell Markers
<b>Form:</b>	Liquid

### Antibody Data

---

<b>Reactivity:</b>	Human
<b>Clone:</b>	R87-5Y-7
<b>Clonality:</b>	Monoclonal Antibody
<b>SwissProt ID:</b>	P04626
<b>Immunogen:</b>	A synthesized peptide derived from human Phospho-ErbB2(Y1221 + Y1222)
<b>Gene ID:</b>	2064

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Gene Name:** ERBB2  
**Host Species:** Rabbit  
**Isotype:** IgG  
**Purification:** Affinity Chromatography  
**Conjugated:** Unconjugated  
**Modification:** Phosphorylated  
**Molecular Weight:** Calculated MW: 138 kDa, Observed MW: 185kDa

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
 Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

	<b>Application</b>	<b>Antibody Dilution Ratio</b>
<b>Antibody Dilution Ratio:</b>	WB	1:500-1:1000

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.