

Anti-MonoMethyl-Histone H3 (Lys14) [R03-7N4] Monoclonal Antibody

AGMB05352

Description

This Anti-MonoMethyl-Histone H3 (Lys14) [R03-7N4] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05352
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	H3K14me, H3.4, H3/g, H3FT, H3t
Applications:	WB ICC/IF
Research Area:	Epigenetics
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse
Clone:	R03-7N4
Clonality:	Monoclonal Antibody
SwissProt ID:	Q16695
Immunogen:	A synthetic peptide of human MonoMethyl-Histone H3-K14
Gene ID:	8290
Gene Name:	H3-4

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Methylated
Molecular Weight:	Calculated MW: 16 kDa, Observed MW: 16 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:2000
	IF	1:50-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.