

Anti-DiMethyl-Histone H3 (Lys23) [5R9-V7-W3] Monoclonal Antibody

AGMB05363

Description

This Anti-DiMethyl-Histone H3 (Lys23) [5R9-V7-W3] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05363
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	H3K23me2, Histone H3/b ; Histone H3/c ; Histone H3/d ; Histone H3/f
Applications:	WB IHC-P
Research Area:	Epigenetics
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse
Clone:	5R9-V7-W3
Clonality:	Monoclonal Antibody
SwissProt ID:	P68431
Immunogen:	Peptide
Gene ID:	8350
Gene Name:	H3C1

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Host Species: Mouse

Isotype: IgG

Purification: Protein G

Conjugated: Unconjugated

Modification: Dimethylated

Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS, Glycerol and BSA

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	IHC-P	1:50-1:100
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.	